

Antibodies Product Datasheet

Product Name:	Acetyl-EPAS-1 (K385) Antibody
Catalog Number:	D12118
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized acetyl-peptide of Acetyl-EPAS-1 (K385)
Concentration:	1 mg/ml
Gene Name:	EPAS1 BHLHE73 HIF2A MOP2 PASD2
Protein Name:	Acetyl-EPAS-1 (K385)
Other Name:	endothelial PAS domain protein 1
Human Swiss Prot No.:	Q99814
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-EPAS-1 (K385) Polyclonal Antibody detects endogenous levels of Acetyl-EPAS-1 (K385)
Storage and Stability:	-20°C/1 year