

Antibodies Product Datasheet

Product Name:	AR(Acetyl-K631) Antibody
Catalog Number:	D12046
Reactivity:	Human:K631, Mouse:K610, Rat:K613
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized Acetyl peptide derived from human AR. at AA range: K631
Concentration:	1 mg/ml
Gene Name:	AR DHTR NR3C4
Protein Name:	AR
Other Name:	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4)
Human Swiss Prot No.:	P10275
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AR(Acetyl-K631), It doesn't react with total protein.
Storage and Stability:	-20°C/1 year