

Antibodies Product Datasheet

Product Name:	NGAL Nano antibody
Catalog Number:	NA0009
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human NGAL.
Concentration:	
Gene Name:	ngal
Protein Name:	
Other Name:	24p3; 25 kDa alpha-2-microglobulin-related subunit of MMP-9; HNL; Lcn2; Lipocalin 2; Lipocalin-2; Migration stimulating factor inhibitor; MSFI; Neutrophil gelatinase-associated lipocalin; NGAL; NGAL_HUMAN; Oncogene 24p3; p25; Siderocalin.
Human Swiss Prot No.:	P80188
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of NGAL.
Storage and Stability:	-20°C/1 year