

Antibodies Product Datasheet

Product Name:	ApoB (Apolipoprotein B Nano antibody
Catalog Number:	NA0008
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human ApoB
Concentration:	
Gene Name:	APOB
Protein Name:	Apolipoprotein B-100
Other Name:	APOB; Apolipoprotein B-100; Apo B-100
Human Swiss Prot No.:	P04114
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of ApoB
Storage and Stability:	-20°C/1 year