
Antibodies Product Datasheet

Product Name:	ApoA-I (Apolipoprotein A-I) Nano antibody
Catalog Number:	NA0007
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human ApoA-I
Concentration:	
Gene Name:	APOA1
Protein Name:	Apolipoprotein A-I
Other Name:	APOA1; Apolipoprotein A-I; Apo-AI; ApoA-I; Apolipoprotein A1
Human Swiss Prot No.:	P02647
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of ApoA-I
Storage and Stability:	-20°C/1 year