

Antibodies Product Datasheet

Product Name:	PA (Prealbumin Nano antibody)
Catalog Number:	NA0006
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human PA
Concentration:	
Gene Name:	TTR PALB
Protein Name:	Prealbumin
Other Name:	Transthyretin (ATTR) (Prealbumin) (TBPA)
Human Swiss Prot No.:	P02766
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of PA
Storage and Stability:	-20°C/1 year