

Antibodies Product Datasheet

Product Name:	PCT (Procalcitonin) Nano antibody
Catalog Number:	NA0005
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human PCT
Concentration:	
Gene Name:	CALCA CALC1
Protein Name:	Calcitonin [Cleaved into: Calcitonin; Katalcalcin (Calcitonin carboxyl-terminal peptide) (CCP) (PDN-21)]
Other Name:	Calcitonin [Cleaved into: Calcitonin; Katalcalcin (Calcitonin carboxyl-terminal peptide) (CCP) (PDN-21)]
Human Swiss Prot No.:	P01258
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of PCT
Storage and Stability:	-20°C/1 year