
Antibodies Product Datasheet

Product Name:	Cystatin c Nano antibody
Catalog Number:	NA0004
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human Cystatin c
Concentration:	
Gene Name:	CST3
Protein Name:	Cystatin-C
Other Name:	CST3; Cystatin-C; Cystatin-3; Gamma-trace; Neuroendocrine basic polypeptide; Post-gamma-globulin
Human Swiss Prot No.:	P01034
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of Cystatin c
Storage and Stability:	-20°C/1 year