
Antibodies Product Datasheet

Product Name:	CD38 Nano antibody
Catalog Number:	NA0003
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human CD38
Concentration:	
Gene Name:	CD38
Protein Name:	ADP-ribosyl cyclase 1
Other Name:	CD38; ADP-ribosyl cyclase 1; Cyclic ADP-ribose hydrolase 1; cADPr hydrolase 1; T10; CD antigen CD38
Human Swiss Prot No.:	P28907
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of CD38
Storage and Stability:	-20°C/1 year