

## Antibodies Product Datasheet

<b>Product Name:</b>	PD-L1 Nano antibody
<b>Catalog Number:</b>	NA0002
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	FCM, ELISA
<b>Dilutions:</b>	ELISA 1:10000-20000
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Nano-antibody with his-tag, Camel
<b>Purification:</b>	The antibody was affinity-purified nickel column using his-tag in antibody
<b>Immunogen:</b>	Recombinant protein of full length human PD-L1
<b>Concentration:</b>	
<b>Gene Name:</b>	CD274 B7H1 PDCD1L1 PDCD1LG1 PDL1
<b>Protein Name:</b>	Programmed cell death 1 ligand 1 (PD-L1) (PDCD1 ligand 1) (Programmed death ligand 1) (B7 homolog 1) (B7-H1) (CD antigen CD274)
<b>Other Name:</b>	Programmed cell death 1 ligand 1 (PD-L1) (PDCD1 ligand 1) (Programmed death ligand 1) (B7 homolog 1) (B7-H1) (CD antigen CD274)
<b>Human Swiss Prot No.:</b>	Q9NZQ7
<b>Formulation:</b>	Liquid in PBS
<b>Specificity:</b>	This antibody detects endogenous levels of PD-L1
<b>Storage and Stability:</b>	-20°C/1 year