
Antibodies Product Datasheet

Product Name:	Anti-COVID-19 Spike mAb IgG (S05FC)
Catalog Number:	CA0004
Reactivity:	COVID-19
Applications:	ELISA ECL Immunogold
Dilutions:	ELISA 1:5000-1:100000, WB 1:1000-1:10000
Modification:	Not Modification Specific
Source:	Humanized mAb IgG
Purification:	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Immunogen:	The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein S
Concentration:	
Gene Name:	SARS-CoV-2 S, 2019-nCoV S, COVID-19 S
Protein Name:	SARS-CoV-2 Spike, 2019-nCoV Spike, COVID-19 Spike
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS
Specificity:	This antibody detects SARS-CoV-2 (2019-nCoV) Spike Protein S1 doman
Storage and Stability:	-20°C/1 year