
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Anti-COVID-19 N protein mAb IgG (S06FC) |
| Catalog Number: | CA0001 |
| Reactivity: | COVID-19 |
| Applications: | ELISA ECL Immunogold |
| Dilutions: | ELISA 1:5000-1:100000, WB 1:1000-1:10000 |
| Modification: | Not Modification Specific |
| Source: | Humanized mAb IgG |
| Purification: | The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen. |
| Immunogen: | The recombinant SARS-CoV-2 (2019-nCoV) N Protein |
| Concentration: | |
| Gene Name: | SARS-CoV-2 N, 2019-nCoV N, COVID-19 N |
| Protein Name: | SARS-CoV-2 Nucleoprotein , 2019-nCoV Nucleoprotein, COVID-19 Nucleoprotein |
| Other Name: | |
| Human Swiss Prot No.: | |
| Formulation: | PBS |
| Specificity: | This antibody detects SARS-CoV-2 (2019-nCoV) Nucleoprotein |
| Storage and Stability: | -20°C/1 year |