

Antibodies Product Datasheet

Product Name:	GDI-1 Antibody
Catalog Number:	R12-2163
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GDI-1. AA range:394-443
Concentration:	1 mg/ml
Gene Name:	GDI1
Protein Name:	Rab GDP dissociation inhibitor alpha
Other Name:	GDI1; GDIL; OPHN2; RABGDIA; XAP4; Rab GDP dissociation inhibitor alpha; Rab GDI alpha; Guanosine diphosphate dissociation inhibitor 1; GDI-1; Oligophrenin-2; Protein XAP-4
Human Swiss Prot No.:	P31150
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GDI-1 Polyclonal Antibody detects endogenous levels of GDI-1 protein.
Storage and Stability:	-20°C/1 year