
Antibodies Product Datasheet

Product Name:	TPPP Antibody
Catalog Number:	C46824
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	TPPP TPPP1
Protein Name:	Tubulin polymerization-promoting protein (TPPP) (25 kDa brain-specific protein) (TPPP/p25) (p24) (p25-alpha)
Other Name:	
Human Swiss Prot No.:	O94811
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	TPPP Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year