
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Rabbit Anti Goat IgG (H&L) - Alexa Fluor 680 |
| Catalog Number: | SA1016 |
| Reactivity: | Goat |
| Conjugate: | Alexa Fluor 680 |
| Applications: | IF/ICC,FCM,WB |
| Dilutions: | Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended. |
| Modification: | Not Modification Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Concentration: | 1mg/mL |
| Other Name: | Rabbit Anti-Goat IgG (H&L) - Alexa Fluor 680 |
| Formulation: | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |
| Storage and Stability: | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |