
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Rabbit Anti Donkey IgG (H&L) - Alexa Fluor 594 |
| Catalog Number: | SA0842 |
| Reactivity: | Donkey |
| Conjugate: | Alexa Fluor 594 |
| Applications: | ELISA,IF/ICC,FCM |
| Dilutions: | IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration) |
| Modification: | Not Modification Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Concentration: | 1mg/mL |
| Other Name: | Rabbit Anti-Donkey IgG (H&L) - Alexa Fluor 594 |
| Formulation: | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |
| Storage and Stability: | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |