
Antibodies Product Datasheet

Product Name:	Goat Anti Human IgG gamma chain - Alexa Fluor 594
Catalog Number:	SA0831
Reactivity:	Human
Conjugate:	Alexa Fluor 594
Applications:	ELISA,IF/ICC,FCM
Dilutions:	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Goat Anti-Human IgG gamma chain - Alexa Fluor 594
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing