

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Dylight 800, Goat Anti Rabbit IgG
<b>Catalog Number:</b>	SA0439
<b>Reactivity:</b>	Rabbit
<b>Conjugate:</b>	Dylight 800
<b>Applications:</b>	IF/ICC,FCM,WB
<b>Dilutions:</b>	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Goat
<b>Purification:</b>	Affinity purified using solid phase goat IgG (H&L) with finally > 95% purity based on SDS-PAGE
<b>Concentration:</b>	
<b>Other Name:</b>	Dylight 800, Goat Anti-Rabbit IgG
<b>Formulation:</b>	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing