
Antibodies Product Datasheet

Product Name:	Dylight 680, Goat Anti Mouse IgG
Catalog Number:	SA0436
Reactivity:	Mouse
Conjugate:	Dylight 680
Applications:	IF/ICC,FCM,WB
Dilutions:	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
Modification:	Not Modification Specific
Source:	Goat
Purification:	Affinity purified using solid phase Mouse IgG (H&L) with finally > 95% purity based on SDS-PAGE
Concentration:	
Other Name:	Dylight 680, Goat Anti-Mouse IgG
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing