
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Dylight 594, Goat Anti Rabbit IgG |
| Catalog Number: | SA0429 |
| Reactivity: | Rabbit |
| Conjugate: | Dylight 594 |
| Applications: | IF/ICC,FCM |
| Dilutions: | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting 1:50-1:1000 dilutions for most fluorescent applications. |
| Modification: | Not Modification Specific |
| Source: | Goat |
| Purification: | Affinity purified using solid phase goat IgG (H&L) with finally > 95% purity based on SDS-PAGE |
| Concentration: | |
| Other Name: | Dylight 594, Goat Anti-Rabbit IgG |
| Formulation: | Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol. |
| Storage and Stability: | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing |