
Antibodies Product Datasheet

Product Name:	DyLight 488, Goat Anti Mouse IgG
Catalog Number:	SA0403
Reactivity:	Mouse
Conjugate:	Dylight 488
Applications:	IF/ICC,FCM
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting 1:50-1:1000 dilutions for most fluorescent applications.
Modification:	Not Modification Specific
Source:	Goat
Purification:	Affinity purified using solid phase Mouse IgG (H&L) with finally > 95% purity based on SDS-PAGE
Concentration:	
Other Name:	DyLight 488, Goat Anti-Mouse IgG
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing