

---

## Antibodies Product Datasheet

<b>Product Name:</b>	DyLight 405, Goat Anti Mouse IgG
<b>Catalog Number:</b>	SA0392
<b>Reactivity:</b>	Mouse
<b>Conjugate:</b>	Dylight 405
<b>Applications:</b>	IF/ICC,FCM
<b>Dilutions:</b>	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting 1:50-1:1000 dilutions for most fluorescent applications.
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Goat
<b>Purification:</b>	Affinity purified using solid phase Mouse IgG (H&L) with finally > 95% purity based on SDS-PAGE
<b>Concentration:</b>	
<b>Other Name:</b>	DyLight 405, Goat Anti-Mouse IgG
<b>Formulation:</b>	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% BSA as stabilizer and 50% Glycerol.
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing