
Antibodies Product Datasheet

Product Name:	Monoclonal mouse Anti Aflatoxin B1
Catalog Number:	SA0021
Reactivity:	AflatoxinB1
Conjugate:	None
Applications:	WB,IHC-p,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:5000-10000 IF 1:5000-10000
Modification:	Not Modification Specific
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Concentration:	4-5mg/mL
Other Name:	Monoclonal mouse Anti-Aflatoxin B1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year