

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Goat anti Rabbit/Mouse polymer HRP, Ready to Use
<b>Catalog Number:</b>	SA0009
<b>Reactivity:</b>	Mouse Rabbit
<b>Conjugate:</b>	HRP
<b>Applications:</b>	IHC-p,IF(paraffin section)
<b>Dilutions:</b>	Ready to Use
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	
<b>Concentration:</b>	
<b>Other Name:</b>	
<b>Formulation:</b>	
<b>Storage and Stability:</b>	4°C/1 year