
Antibodies Product Datasheet

Product Name:	BIOTIN* Goat Anti Mouse IgG(H+L)
Catalog Number:	SA0005
Reactivity:	Mouse
Conjugate:	Biotin
Applications:	WB,IHC,IF/ICC,ELISA
Dilutions:	WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:500-1000
Modification:	Not Modification Specific
Source:	Goat
Purification:	The antibody was affinity-purified from goat antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration:	4 mg/ml
Other Name:	BIOTIN* Goat Anti-Mouse IgG(H+L)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage and Stability:	-20°C/1 year