

Antibodies Product Datasheet

Product Name:	PDGFR beta (Phospho-Tyr740) Antibody
Catalog Number:	P12-1162
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	PDGFR beta (Phospho-Tyr740) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human PDGFR beta around the phosphorylation site of Tyr740.
Concentration:	1 mg/ml
Gene Name:	PDGFRB
Other Name:	CD140b; kinase PDGFR-beta; PDGF-R-beta; PDGFR; PDGFR-beta; PGFRB; platelet-derived growth factor receptor, beta polypeptide
Human Swiss Prot No.:	P09619
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PDGFR beta (Phospho-Tyr740) Antibody detects endogenous levels of PDGFRB only when phosphorylated at Tyr740.
Storage and Stability:	Stable at -20°C for at least 1 year.