

Antibodies Product Datasheet

Product Name:	SREBP-1 (Phospho-Ser439) Antibody
Catalog Number:	P12-1181
Reactivity:	Human, Rat
Applications:	WB, IHC, ELISA
Dilutions:	SREBP-1 (Phospho-Ser439) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human SREBP-1 around the phosphorylation site of Ser439.
Concentration:	1 mg/ml
Gene Name:	SREBF1
Other Name:	SRBP1; SREBF1; SREBP1; sterol regulatory element-binding protein 1; sterol regulatory element-binding transcription factor 1
Human Swiss Prot No.:	P36956
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	SREBP-1 (Phospho-Ser439) Antibody detects endogenous levels of SREBF1 only when phosphorylated at Ser439.
Storage and Stability:	Stable at -20°C for at least 1 year.