

Antibodies Product Datasheet

Product Name:	Survivin (Phospho-Thr34) Antibody
Catalog Number:	P12-1183
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Survivin (Phospho-Thr34) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human Survivin around the phosphorylation site of Thr34.
Concentration:	1 mg/ml
Gene Name:	BIRC5
Other Name:	API4; apoptosis inhibitor 4; apoptosis inhibitor survivin; baculoviral IAP repeat-containing protein 5; BIRC5; IAP4; TIAP
Human Swiss Prot No.:	O15392
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Survivin (Phospho-Thr34) Antibody detects endogenous levels of BIRC5 only when phosphorylated at Thr34.
Storage and Stability:	Stable at -20°C for at least 1 year.