

Antibodies Product Datasheet

Product Name:	TOP2A (Phospho-Ser1525) Antibody
Catalog Number:	P12-1191
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	TOP2A (Phospho-Ser1525) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human TOP2A around the phosphorylation site of Ser1525.
Concentration:	1 mg/ml
Gene Name:	TOP2A
Other Name:	alpha isozyme; DNA topoisomerase II; DNA topoisomerase II-alpha; P11388-1; TOP-2; TOP2; TOP2A; TP2A
Human Swiss Prot No.:	P11388
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TOP2A (Phospho-Ser1525) Antibody detects endogenous levels of TOP2A only when phosphorylated at Ser1525.
Storage and Stability:	Stable at -20°C for at least 1 year.