
Antibodies Product Datasheet

Product Name:	p130 Cas Antibody
Catalog Number:	R12-2284
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	p130 Cas Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p130 Cas.
Concentration:	1 mg/ml
Gene Name:	BCAR1
Other Name:	Breast cancer anti-estrogen resistance protein 1; CRK-associated substrate; Cas scaffolding protein family member 1; p130cas; BCAR1; CAS; CASS1; CRKAS
Human Swiss Prot No.:	P56945
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	p130 Cas antibody detects endogenous levels of p130 Cas.
Storage and Stability:	Stable at -20°C for at least 1 year.