

Antibodies Product Datasheet

Product Name:	PKD1/PKC-mu Antibody
Catalog Number:	R12-2311
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	PKD1/PKC-mu Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PKD1.
Concentration:	1 mg/ml
Gene Name:	KPCD1
Other Name:	KPCD1; PKC-mu; PKCM; PKD; PRKCM; PRKD1; Protein kinase C, mu type; Protein kinase D; kinase PKD1; nPKC-mu
Human Swiss Prot No.:	Q15139
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PKD1/PKC-mu antibody detects endogenous levels of PKD1/PKC-mu.
Storage and Stability:	Stable at -20°C for at least 1 year.