

Antibodies Product Datasheet

Product Name:	ALDOA Antibody
Catalog Number:	R12-2476
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	ALDOA Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ALDOA.
Concentration:	1 mg/ml
Gene Name:	ALDOA
Other Name:	ALDA; ALDOA; ALFA; Fructose-bisphosphate aldolase A; Lung cancer antigen NY-LU-1; Muscle-type aldolase
Human Swiss Prot No.:	P04075
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ALDOA Antibody antibody detects endogenous levels of ALDOA.
Storage and Stability:	Stable at -20°C for at least 1 year.