

Antibodies Product Datasheet

Product Name:	B-RAF Antibody
Catalog Number:	R12-2545
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	B-RAF Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human B-RAF.
Concentration:	1 mg/ml
Gene Name:	BRAF1
Other Name:	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; B-Raf proto-oncogene serine/threonine-protein kinase; BRAF; BRAF1; c-RMIL; kinase B-Raf; p94; RAFB1; RMIL; RMIL serine/threonine-protein kinase; v-Raf murine sarcoma viral oncogene homolog B1
Human Swiss Prot No.:	P15056
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	B-RAF Antibody antibody detects endogenous levels of B-RAF.
Storage and Stability:	Stable at -20°C for at least 1 year.