

Antibodies Product Datasheet

Product Name:	MADD Antibody
Catalog Number:	R12-2984
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	MADD Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MADD.
Concentration:	1 mg/ml
Gene Name:	MADD
Other Name:	MAP kinase-activating death domain protein; Differentially expressed in normal and neoplastic cells; Insulinoma glucagonoma clone 20; Rab3 GDP/GTP exchange factor; DENN; IG20; KIAA0358
Human Swiss Prot No.:	Q8WXG6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	MADD Antibody antibody detects endogenous levels of MADD.
Storage and Stability:	Stable at -20°C for at least 1 year.