

Antibodies Product Datasheet

Product Name:	TGF beta Receptor I Antibody
Catalog Number:	R12-3618
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	TGF beta Receptor I Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TGF beta Receptor I.
Concentration:	1 mg/ml
Gene Name:	TGFBR1
Other Name:	ADRB2; ALK-5; Activin receptor-like kinase 5; ESK2; SKR4; TGF-beta receptor type I precursor; TGF-beta type I receptor; TGF-betaR1; TGFR-1; TGFR1; kinase TGF-beta-R1; receptor R4
Human Swiss Prot No.:	P36897
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TGF beta Receptor I Antibody antibody detects endogenous levels of TGF beta Receptor I.
Storage and Stability:	Stable at -20°C for at least 1 year.