
Antibodies Product Datasheet

Product Name:	ZDHHC7 Antibody
Catalog Number:	R12-3691
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	ZDHHC7 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZDHHC7.
Concentration:	1 mg/ml
Gene Name:	ZDHHC7
Other Name:	Palmitoyltransferase ZDHHC7; Zinc finger DHHC domain-containing protein 7; DHHC-7; Zinc finger protein 370
Human Swiss Prot No.:	Q9NXF8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZDHHC7 Antibodyantibody detects endogenous levels of ZDHHC7.
Storage and Stability:	Stable at -20°C for at least 1 year.