

Antibodies Product Datasheet

Product Name:	ZFYVE20 Antibody
Catalog Number:	R12-3693
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	ZFYVE20 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZFYVE20.
Concentration:	1 mg/ml
Gene Name:	ZFYVE20
Other Name:	Rabenosyn-5; FYVE finger-containing Rab5 effector protein rabenosyn-5; Zinc finger FYVE domain-containing protein 20; 110 kDa protein Q9H1K0
Human Swiss Prot No.:	Q9H1K0
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZFYVE20 Antibody antibody detects endogenous levels of ZFYVE20.
Storage and Stability:	Stable at -20°C for at least 1 year.