
Antibodies Product Datasheet

Product Name:	ZNF23 Antibody
Catalog Number:	R12-3695
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	ZNF23 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF23.
Concentration:	1 mg/ml
Gene Name:	ZNF23
Other Name:	Zinc finger protein 23; Zinc finger protein 359; Zinc finger protein KOX16; ZNF23; KOX16; ZNF359
Human Swiss Prot No.:	P17027
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZNF23 Antibody antibody detects endogenous levels of ZNF23.
Storage and Stability:	Stable at -20°C for at least 1 year.