
Antibodies Product Datasheet

Product Name:	ZNF287 Antibody
Catalog Number:	R12-3697
Reactivity:	Human
Applications:	WB, IHC, ELISA
Dilutions:	ZNF287 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZNF287.
Concentration:	1 mg/ml
Gene Name:	ZNF287
Other Name:	Zinc finger protein with KRAB and SCAN domains 13
Human Swiss Prot No.:	Q9HBT7
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZNF287 Antibody antibody detects endogenous levels of ZNF287.
Storage and Stability:	Stable at -20°C for at least 1 year.