

Antibodies Product Datasheet

Product Name:	Anti-Human Androgen Receptor (ABT101) Mouse Monoclonal Antibody
Catalog Number:	V0179
Clone ID:	ABT101
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	AR
Gene Name:	P10275
Other Name:	AIS, Androgen receptor, AR, AR8, DHTR, Dihydro testosterone receptor, Dihydrotestosterone receptor (DHTR), HUMARA, HYPSP1, KD, NR3C4, Nuclear receptor subfamily 3 group C member 4 (NR3C4), Nuclear receptor subfamily 3 group C member 4, SBMA, SMAX1, Spinal and bulbar muscular atrophy (SBMA), Spinal and bulbar muscular atrophy, Testicular Feminization (TFM), TFM
Isotype	IgG2b, kappa
Localization:	Prostate
Positive Tissue Control:	Nuclear
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human AR protein.
Storage and Stability:	-20°C/1 year