

Antibodies Product Datasheet

Product Name:	CDX-2 (PT0062) mouse mAb
Catalog Number:	N1827
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG3, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CDX-2
Concentration:	
Gene Name:	CDX2 CDX3
Protein Name:	Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2)
Other Name:	
Human Swiss Prot No.:	Q99626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CDX-2. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year