

## Antibodies Product Datasheet

<b>Product Name:</b>	Oct-2 (PT0202) mouse mAb
<b>Catalog Number:</b>	N1770
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	IHC-p,IF(paraffin section)
<b>Dilutions:</b>	IHC-p 1:100-500
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse/IgG1, Kappa
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human Oct-2
<b>Concentration:</b>	
<b>Gene Name:</b>	POU2F2 OCT2 OTF2
<b>Protein Name:</b>	POU domain, class 2, transcription factor 2 (Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2) (Octamer-binding protein 2) (Oct-2) (Octamer-binding transcription factor 2) (OTF-2)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P09086
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of human Oct-2. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
<b>Storage and Stability:</b>	-20°C/1 year