
Antibodies Product Datasheet

Product Name:	Estrogen Receptor β mouse mAb(PT0318)
Catalog Number:	N1747
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),WB
Dilutions:	IHC-p 1:100-300, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Estrogen Receptor β
Concentration:	1 mg/ml
Gene Name:	ESR2 ESTRB NR3A2
Protein Name:	Estrogen Receptor β
Other Name:	Estrogen receptor beta (ER-beta) (Nuclear receptor subfamily 3 group A member 2)
Human Swiss Prot No.:	Q92731
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Estrogen Receptor β
Storage and Stability:	-20°C/1 year