
Antibodies Product Datasheet

Product Name:	Gamma Tubulin mouse Monoclonal Antibody(7B1)
Catalog Number:	N1634
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Gamma Tubulin
Concentration:	1 mg/ml
Gene Name:	TUBG1
Protein Name:	TUBG1
Other Name:	TUBG1
Human Swiss Prot No.:	P23258
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gamma Tubulin protein detects endogenous levels of TUBG1
Storage and Stability:	-20°C/1 year