
Antibodies Product Datasheet

Product Name:	Muscle actin mouse Monoclonal Antibody(4G10)
Catalog Number:	N1626
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Muscle actin
Concentration:	1 mg/ml
Gene Name:	ACTA1
Protein Name:	ACTA1
Other Name:	ACTA1
Human Swiss Prot No.:	P68133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Muscle actin protein detects endogenous levels of ACTA1
Storage and Stability:	-20°C/1 year