

Antibodies Product Datasheet

Product Name:	CREB-1 Monoclonal Antibody(4B4)
Catalog Number:	N1516
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF (paraffin section)
Dilutions:	WB: 1:500-1000 IHC: 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of CREB-1
Concentration:	
Gene Name:	CREB1
Protein Name:	Cyclic AMP-responsive element-binding protein 1
Other Name:	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1
Human Swiss Prot No.:	P16220
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CREB-1 protein.
Storage and Stability:	-20°C/1 year