

Antibodies Product Datasheet

Product Name:	CDX-2 Mouse mAb(8D9)
Catalog Number:	N1501
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CDX-2
Concentration:	1 mg/ml
Gene Name:	CDX2 CDX3
Protein Name:	CDX-2
Other Name:	Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2)
Human Swiss Prot No.:	Q99626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.74% sodium azide.
Specificity:	This antibody detects endogenous levels of CDX-2 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year