

Antibodies Product Datasheet

Product Name:	EGFR Monoclonal Antibody(M6)
Catalog Number:	N1399
Reactivity:	Human
Applications:	WB
Dilutions:	WB: 1:1000-3000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of EGFR
Concentration:	
Gene Name:	EGFR
Protein Name:	Epidermal growth factor receptor
Other Name:	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
Human Swiss Prot No.:	P00533
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous EGFR protein.
Storage and Stability:	-20°C/1 year