

Antibodies Product Datasheet

Product Name:	CDX2 Monoclonal Antibody(14H6)
Catalog Number:	N1322
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, IHC-p
Dilutions:	WB: 1:1000 IHC: 1:200 IF 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CDX2
Concentration:	
Gene Name:	CDX2
Protein Name:	Homeobox protein CDX-2
Other Name:	CDX2; CDX3; Homeobox protein CDX-2; CDX-3; Caudal-type homeobox protein 2
Human Swiss Prot No.:	Q99626
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CDX2 proteins.
Storage and Stability:	-20°C/1 year