

## Antibodies Product Datasheet

<b>Product Name:</b>	EGFR Monoclonal Antibody(1B10)
<b>Catalog Number:</b>	N1308
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IHC-p,IF/ICC
<b>Dilutions:</b>	WB: 1:1000 IHC-p: 50-200 IF: 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of EGFR
<b>Concentration:</b>	
<b>Gene Name:</b>	EGFR
<b>Protein Name:</b>	Epidermal growth factor receptor
<b>Other Name:</b>	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
<b>Human Swiss Prot No.:</b>	P00533
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous EGFR proteins.
<b>Storage and Stability:</b>	-20°C/1 year